# geko™

T-2

circulation support

## EN

### Warnings

Implanted electronic devices (for example a cardiac pacemaker): - specialist medical opinion should be obtained before the patient is subjected to stimulation. Recently diagnosed or suspected DVT; specialist medical opinion should be obtained before the patient is subjected to stimulation.

Read the instructions carefully before use. Only apply the device in accordance with instructions. Improper use of this device can be dangerous.

The geko™ device may be operated either by the patient or a third party. Patients should be trained in the positioning and operation of the device.

Keep out of the reach of children and pets.

No modification of this equipment is allowed.

Do not use whilst driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user or others at undue risk of injury.

Stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.

Application of electrodes near the thorax may increase the risk of cardiac fibrillation.

Do not apply over sore, infected or inflamed areas, broken skin or skin eruptions, e.g. phlebitis, thrombophlebitis, varicose veins etc.

Do not apply over or, in proximity to, cancerous lesions.

Do not use in proximity (i.e. within 1m) of short wave/ microwave equipment, as this may affect the device.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Do not use if packaging is open or appears damaged.

Do not use a device that appears to be damaged.

Do not use the device in close proximity to heat sources, such as fires or radiant heaters, as excess heat may affect the performance of the device.

In some cases skin inflammation or irritation can develop in the contact area: either remove the device or re-attach in the alternative fitting location. If the condition persists or recurs obtain specialist medical advice before resuming use.

Medical advice should be sought when the patient is pregnant, has diagnosed heart conditions or epilepsy, or following surgery where muscle contractions may disrupt the healing process.

The device should be removed before the patient

undergoes MRI as it contains ferromagnetic components.

The device should be switched off during ECG monitoring using leg electrodes as it may interfere with ECG leg electrode signals.

The device is not intended for use in an oxygen rich environment.

### Precautions

Check the device is functioning correctly every 2 hours.

Use the device only as instructed. Visit [www.gekodevices.com](http://www.gekodevices.com) for further advice and usage tips. Keep out of the reach of children.

Do not switch on unless correctly attached to the patient.

The device must be kept clean and dry. Do not bath or shower whilst wearing the device.

Do not apply a device that has been worn by another individual as this will carry a risk of infection.

Specialist medical opinion should be obtained before devices are used for more than 28 consecutive days.

Remove the device carefully, lifting the long-end first and applying a constant gentle force.

Excessive force may damage the device.

### EMC Declaration

The geko™ device is certified to EN 60601-1-2:2007 regarding Electromagnetic Compatibility. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided at: <http://www.gekodevices.com/>

### Disposal

Batteries must be disposed of in accordance with any local legislation. The battery can be removed by breaking open the housing and prising out.

See [www.gekodevices.com](http://www.gekodevices.com) for guidance.

Do not incinerate the device.

### Symbols

⊕ On/Increase stimulation button ⊖ Off/decrease stimulation button

### Firstkind Limited\*

Hawk House
Peregrine Business Park
High Wycombe
HP13 7DL
United Kingdom
+44 (0)1494 572040

### www.gekodevices.com

\*Firstkind Limited is a wholly owned subsidiary of Sky Medical Technology Limited.

If assistance is required or to report difficulties, unexpected operation or events, contact the manufacturer.

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**geko™** and **OnPulse™** are registered trademarks of Sky Medical Technology Limited. This product is covered by limited and pending patents. ©2017 Sky Medical Technology Limited.
\*Firstkind Limited er et helejet datterselskab af Sky Medical Technology Limited.